



# THE DETERMINATION OF LACTIC ACID IN MILK :

## COMPARISON BETWEEN FRESH MILK AND MILK LEFT IN THE OPEN AIR FOR 3 DAYS

**Objective :** to compare fresh milk and fermented milk by a simple determination of the lactic acid  $CH_2-CH(OH)-COOH$  with a soda solution .

**Note :** the acidity of milk is conventionally expressed in Domic degrees : one Domic degree °D corresponds to 0,1g of lactic acid per litre of milk .



### Materials and products :

- fresh milk and fermented milk
- distilled water, phenolphthalein soda solution of concentration  $C_0=0.05M$
- a magnetic stirrer, 2 Erlenmeyer flasks, a burette, volumetric pipettes (10mL, 20mL), pH meter, volumetric flask (50 mL), a test tube, a funnel and a stand

### Procedure :

- take 20mL of milk in a volumetric pipette,
- introduce this volume into a 200mL Erlenmeyer flask, then complete with 100mL of distilled water,
- measure the pH,
- add a few drops of phenolphthalein,
- add to the burette a soda solution  $C_0=0.05M$ . Let  $V_{e1}$  be the volume poured at the equivalence.



$$M(\text{lactic acid}) = \frac{C_0 \times V_{eq} \times 89}{20}$$